

Product datasheet for **RR203245**

Raet1l (NM_001013063) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Raet1l (NM_001013063) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Raet1l
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR203245 representing NM_001013063 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCAAAGAGAGGAACCAAGTGAATCTTTCTCTGAGTCTGGCCCTTTGCTGCTGCTGAATTTTC
TGGGACCGCGATGCTAGAAGATGTACACTCTCTTAGTTGCACTTGATCATCAGGTCTCCTACCAAAGC
AGAACCTCTCTGGTATGAAGGACAGTGTCAATGGATGGAATTGATTTCTTTCTTTAATAATGTGAAC
AAGACCATGCCTTTGGGTAAACCAGGAGAAGATGGCAAATGCCACCGAAGTGTGGACATGTGTGCCCAAGT
GGATAGACCTGTTGGGCCAGGACTTAAGGGACAAGGTGTCTGAAATCAAAGGGGACACTACCAAGACCAG
TGGTCACCCCACTTTATGGGTCAACATGCTTTCTAACCACAGCCATAGACATATCATTGGTGCTTCTCTGT
GTATTC AACATCAGTGAGAATTATTCCTTCATCTTGGACCTAAAGACTATGAAGTGGAACTGATTAATC
CTGTACCAGAGGACATAATGAAAGAATGAAAGATGATGCGACATTAATGAAACATCTGAATTTGTCCAT
AGCAGAATGCAGTCAGAAATTCATGAATTCTTAAATCAGCACGAGGAAAGAACAAGATCAAGAACAAGG
GCCCCAGATATCGCCAGCCTCCACCTACCAGGCAGTTTCCTTAAGGAAGTATTGATCCCTGTGGTTC
TCATCATCTTTATTTTTCATTTGCATATGTGTGAAGAGGATGTGTTGTACCCAAGGAGGTAGAGA
ACTTCAGCTGCTCAAGGACTACTACATCAACCTCCACAACCAGAGGAATGAGGGGCTAGAGGAGACAAAA
GGAGAGGTTCCCAACCAGAGCACCAAAATGTTGTCCCAGGAGGAAGATCCCAAGACATTGGGAGGCCCAA
GGAGGGCATCAGAATCTGTGGAAGTGGAGTTACAGATGTTTTCAAGTACTGGGAGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR203245 representing NM_001013063
Red=Cloning site Green=Tags(s)

MAKRGTTKCNLSLSLGLLLLLNFLGTAMLEDVHSLSCNLIIRSPTKAEP LWYEGQCSMDGIDFLSFNNVN
 KTMP L GNQEKMANATEVWTCVTQWIDL LGQDLRDKVSEIKGDTTKTSGHPTLWVTMLS NHSHRHIIGASC
 VFNI SENYSFILD LKTMKWKLINPVPEDIMKEWKDDATLMKHLNLSIAECSQKFNEFLNQHEERTSRTR
 APDIAQPPPT RQFPYKEVLI PVVLIIFII FICIVKRMCTQGGRELQLLKDYIYNLHKPQRNEGLEETK
 GEV PNQSTKM L SQEEDPKTLGGPRRASESVELELQMFSS TGS


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001013063

ORF Size:	966 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001013063.1 , NP_001013081.1
RefSeq Size:	1459 bp
RefSeq ORF:	969 bp
Locus ID:	292461
Cytogenetics:	7q34
MW:	36.6 kDa