

## Product datasheet for **RR203250L3V**

### Nr3c1 (NM\_012576) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Nr3c1 (NM_012576) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Nr3c1
Synonyms:	Gcr; GR; Grl
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_012576
ORF Size:	2382 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR203250).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_012576.2</a> , <a href="#">NP_036708.2</a>
RefSeq Size:	6327 bp
RefSeq ORF:	2385 bp
Locus ID:	24413
UniProt ID:	<a href="#">P06536</a>
Cytogenetics:	18p11
Gene Summary:	glucocorticoid receptor that binds and activates hormone-dependent transcriptional enhancers [RGD, Feb 2006]



[View online »](#)