

Product datasheet for **RR203267L3V**

Ccin (NM_001005904) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ccin (NM_001005904) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ccin
Synonyms:	MGC94933
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001005904
ORF Size:	1764 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR203267).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001005904.1 , NP_001005904.1
RefSeq Size:	1919 bp
RefSeq ORF:	1767 bp
Locus ID:	298392
UniProt ID:	Q5XI58
Cytogenetics:	5q22
Gene Summary:	Possible morphogenetic cytoskeletal element in spermiogenic differentiation. [UniProtKB/Swiss-Prot Function]



[View online »](#)