

Product datasheet for **RR203284L4V**

Galnt10 (NM_130742) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Galnt10 (NM_130742) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Galnt10
Synonyms:	Galnt9
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_130742
ORF Size:	1809 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR203284).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_130742.1 , NP_570098.1
RefSeq Size:	1812 bp
RefSeq ORF:	1812 bp
Locus ID:	170501
UniProt ID:	Q925R7
Cytogenetics:	10q22
Gene Summary:	UDP-GalNAc:polypeptide N-acetylgalactosaminyltransferase; catalyzes the transfer of GalNAc from the nucleotide sugar UDP-GalNAc to the hydroxyl group of either serine or threonine [RGD, Feb 2006]



[View online »](#)