

Product datasheet for **RR203319**

Cnot2 (NM_001011988) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cnot2 (NM_001011988) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Cnot2
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR203319 representing NM_001011988 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGGGGGATGAGCAACAATACCCCTCAGTTAAATCGCAGCTTATCACAAGGCACTCAGTTACCGAGCC
ACGTCACGCCAACACAGGGGTACCAACAATGTCACTTCACACGCCTCCATCTCCAAGCAGGGGTATTTT
GCCTATGAATCCTAGGAATATGATGAACCACTCCAGGTTGGTCAGGGCATTGGAATTCCTAGCAGGACA
AATAGCATGAGCAGTTCAGGGTTAGGTAGCCCAACAGAAGCTCGCCAAGCATAATATGATGCCAAAGC
AGCAGCCTTCTCGACAGCCTTTTACTGTGAACAGTATGTCTGGATTTGGAATGAACAGGAATCAGGCATT
TGGAAATGAATAACTCCTTATCAAGTAACATTTTAAATGGAACAGATGGCAGTGAAAATGTGACAGGATTG
GACCTTTCAGATTTCCAGCATTAGCAGACCGGAATAGAAGGGAAGGAAGTGGAACCCAACTCCATTAA
TAAACCCCTTGGCTGGAAGAGCCCTTACGTTGGAATGGTAACAAAACCTGCAAATGAGCAATCCAGGA
CTTCTCAATACACAACGAAGATTTCCAGCACTACCTGGCTCCAGCTATAAAGACCCAACGTCGAGTAAT
GACGACAGCAAGTCTAATTTGAGTACATCGGGGAAGCAGCTTCAAGTACAGATGGGCCAAATTTCTCTG
GAGATAAAAGTTCAACAACACAAAACAATAATCAGCAGAAAAAGGGGATCCAGGTGTACCTGATGGTCCG
AGTTACTAACATTCCTCAAGGGATGGTGACGGACCAATTTGGAATGATTGGCCTGTTAACATTTATCAGG
GCAGCAGAGACAGACCCGGGAATGGTACATCTTGCTTTAGGAAGTGACTTAACAACATTAGGCCTCAATC
TGAACACCCGAAAATCTTACCCAAATTTGCATCACCTGGGCTTCTTACCTTGTGACCTCAAGA
TATAGACTTCCATGTTCCGTGAGAATTTAAACAAACATTCACATTAGGGATAAGCTGGCTGCAATCAA
CTTGGCCGATATGGAGAAGACCTCCTCTTCTATCTTATTACATGAACGGAGGAGACGTATTACAACTCT
TAGCTGCAGTAGAGCTTTTTAACCGTATTGGAGATACCACAAGGAGGAGCGGGTATGGATTACCAGGGC
ACCAGGCATGGAGCCAACGATGAAAACCAACACATATGAGCGAGGAACCTACTACTTCTTTGACTGTCTC
AACTGGAGGAAAGTAGCTAAGGAACAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR203319 representing NM_001011988
 Red=Cloning site Green=Tags(s)


MRGMSNNTPQLNRSLSQGTQLPSHVTPPTGVPTMSLHTPPSPSRGILPMNPRNMMNHSQVGQIGIPSR
 NSMSSSGLGSPNRSSPSIICMPKQQPSRQPFVNSMSGFGMNRNQAFGMNSSLSSNIFNGTDGSENV
 TGLDLSDFPALADRNRREGSGNPTPLINPLAGRAPYVGMVTKPANEQSQDFSIHNEDFPALPGSSYK
 DPTSSNDDSKSNLSTSGKTTSSTDGPKFPGDKSSTTQNNNQKKGIVLQVLPDGRVTNIPQGMVTDQ
 FGMIGLLTFIRAAETDPGMVHLALGSDLTTLGLNLNSPENLYPKFASPWASSPCRQDIDFHPVSEY
 LTNIHIRDKLAAIKLGRYGEDLLFYLYMNGGDVLLA AVELFNRDWRYHKEERVWITRAPGMEPTMK
 TNTYERGTYFFDCLNWRKVAKEH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:	NM_001011988
ORF Size:	1287 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001011988.1 , NP_001011988.1
RefSeq Size:	1810 bp
RefSeq ORF:	1290 bp
Locus ID:	299805
Cytogenetics:	7q22
MW:	47.1 kDa