

Product datasheet for **RR203361L4V**

Tap2 (NM_032056) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tap2 (NM_032056) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Tap2
Synonyms:	Abcb3; Cim
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_032056
ORF Size:	2109 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR203361).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_032056.3 , NP_114445.2
RefSeq Size:	2446 bp
RefSeq ORF:	2112 bp
Locus ID:	24812
UniProt ID:	P36372
Cytogenetics:	20p12
Gene Summary:	transports peptides into the ER lumen for binding with MHC class I molecules; plays a role in antigen processing and presentation [RGD, Feb 2006]



[View online »](#)