

Product datasheet for **RR203374**

Tssc4 (NM_001013194) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tssc4 (NM_001013194) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Tssc4
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR203374 representing NM_001013194
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCAGGTCCTTATGGCGGAGACAGAGGCCGGCTTGGAGGTTGAGGAACCCACCGAGGATGACACCCTGC
CTTCTGACACCGTCTCCCTCAGCGACTCGGACTCTGACCTCAGTTACCTAGTGGTGTGGAGTCCAAGC
GTTGTCCCAGAGAGGCTTCTGGGGAGAGCCAGGAGGACTCGGGCCCTGATGACCTCCCTCGCACCCC
ACGGGCATTCCCACCACTGCAGTGCAGCCGTTCCACCTCCGAGGCATGAGCTCCACTTCTCCCAGCGGA
GCCACAGCATCTTTGATTGTCTGGAGAGTGCAGCCCGGCAGGAACCGTGCTCTGCTCCCCAAACAGTGT
GGTTGACAATTGCAGTCAAGCGGCCTGTGGCTCCCCAAGTCAGACTCCAGCAAGAAGCTTGAGCAGA
GTGCATGGGAACACTGACCCGACCAGGGTACACCCTGTCCCTGACTATGTGTCACACCCTGAGCGTTGGA
CCAAATACAGTCTGGAAGATGTGTCTGAAACCAGTGCAGCAGCAATCGTGATGCTGCCCTGGCCTTCT
GAGCTCTCGAAGCCAGGCATCCCCACAGACTATGTGCCCTTTTCAACCAGGACCCCTAGCTGTGGG
GAGGGAAGAGTGGTCTTACCAAACAGTGCAGGGCAGCGAGGCTAGGGCTGAGAGGAAGAGGGTCTAA
AGAAGGGGTTGTGTGAGGTGCAGGGGTGAGGCCTCAGTGGAGCTAGCCATCTCGCCGGGCTGAGGC
TGAGGAATGGAGTGGCCACCAGGACAGCCAGAGTGGTGTAGTGCCTTCAAGAGCTGCCCTGCTGAGTCC
TCATCAGGCCCATAGGTATGAAAAGTGTGGTTCCATGGCAGCAAGAAGCGGAGCAGAGACCCTTCC
GGAACAGGGATAGTAACCCGGAGGGGCCAGGGTCTGAGAGAGGACCCTCAGTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR203374 representing NM_001013194
 Red=Cloning site Green=Tags(s)

MSGLMAETEAGLEVEEPTEDDTLPSTVSLSDSDSLPSGVEVQALSPERLSGESQEDSGPDDPPSHP
 TGIPTTAVQPFFHLRGMSTFSQRSHSIFDCLESAARQEPESAPQTSVVDNCSFKRPVAPPSQTPARLSR
 VHGNTDPTRVHPVPDYVSHPERWTKYSLEDVSETSEQSNRDAALAFSSRSQASPTDYVPPFNQDPSSCG
 EGRVVF TKPVRGSEARAERKRVLKKGVVSGAGGEASVELAHLAGPEAEWWSGHQGPVEVVPSEARPES
 SSGPIGMKTVGFHGSKKRSRDHFRNRDSNPEGPGSERGPSV


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001013194

ORF Size:	963 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001013194.2 , NP_001013212.2
RefSeq Size:	1451 bp
RefSeq ORF:	966 bp
Locus ID:	361682
Cytogenetics:	1q42
MW:	34.3 kDa