

## Product datasheet for **RR203428**

### Fli1 (NM\_001017381) Rat Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Fli1 (NM_001017381) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Fli1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR203428 representing NM_001017381 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**C

ATGGACGGGACTATTAAGGAGGCTCTGTCAGTGGTGGAGTGACGACCAGTCCCTTTTTGATTACGCATACG  
GAGCGGCCGCCATCTCCCAAGGCAGATAGACTGCTTCGGGGAGTCTGACTACGGGCAGCCCCACAA  
AATCAACCCCTGCCACCACAGCAGGAGTGGATCAACCAGCCAGTGAAGGTCAATGTCAAGCGGAATAT  
GATCACATGAATGGATCCAGGGAGTCTCCGGTGGACTGCAGTGTGAGCAATGCAACAAGCTGGTGGGTG  
GAGGCGAAGCCAACCCATGAACTACAATAGCTACATGGATGAGAAGAATGGCCCCCTCTCCCAACAT  
GACCACCAACGAGCGGAGAGTCAATTGTGCCTGCAGACCCACACTGTGGACACAGGAGCAGTTCGACAG  
TGGCTGGAGTGGCCATAAAGGAATATGGCTTGTGGAGATTGACACTTCCTTCTCCAGAACATGGATG  
GCAAGGAAGTGTGAAAATGAGCAAGGAGGACTTCCTCCGAGCCACCTCCCTCTACAACACAGAAGTGT  
GTTGTGCGACCTCAGTTACCTCAGGAAAGTTCAGTCTGCTGAGGAGGGGGATGGAGCAACAATA  
TGAACTCTGGCATCAACAAAAGTCTCCCTTGGAGGAGCACAGAGCATGGGCAAGAACACTGAGCAACG  
GCCCCAGCCAGATCCTTATCAGATCCTGGGGCCACCAGCAGCCGCTAGCAAACCTGGAAGTGGGCAG  
ATCCAAGTGTGGCAGTTTCTCCTGGAAGTACTGTCTGACAGCGCAACGCCAGTTGCATTACTTGGGAGG  
GGACCAACGGGAATTCAAAATGACGGACCCTGATGAGGTGGCCAGGCGCTGGGGAGAGCGGAAGAGCAA  
GCCCAACATGAATTACGATAAGCTGAGCCGGGCCCTCCGATACTACTACGATAAAAAATTATGACCAAA  
GTGCATGGCAAAAGGTATGCCTACAAGTTTGACTTCCATGGCATTGCCAAAGCCCTGCAGCCACACCCAA  
CAGAGACATCAATGTACAAGTATCCTTCTGACATCTCCTACATGCCTTCTACCATGCCATCAACAAAA  
GGTGAACCTTGTCCAGCCCCACCCATCTCCATGCCTGTACCTCCTCCAGCTTCTTTGGAGCAGCGTCA  
CAATACTGGACCTCCCCACTGCGGGATCTACCCGAACCCAGTGTACCCGCCATCCTAACCCACAG  
TGCCCTCACACTTAGGCAGCTACTAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RR203428 representing NM\_001017381  
Red=Cloning site Green=Tags(s)


MDGTIKEALSVVSDQSLFDSAYGAAHLPKADMTASGSPDYGQPHKINPLPPQEQWINQPVRVNVKREY  
 DHMNGSRESPVDCSVSKCNKL VGGGEANPMNYSYMDEKNGPPPNMTTNERRVI VPADPTLWTQEHVRQ  
 WLEWAIKEYGLMEIDTSFFQNMDGKELCKMSKEDFLRATSLYNTVEVLLSHLSYLRESSLLAYNTPSHTDQ  
 SSRLNVKEDPSYDSVRRGGWSNNMNSGINKSPPLGGAQSMGKNTQRPQPDYQILGPTSSRLANPGSGQ  
 IQLWQFLLELLSDSANASCITWEGTNGEFKMTDPDEVARRWGERKSKPNMNYDKLSRALRYYYDKNIMTK  
 YHGKRYAYKDFDHGIAQALQPHPTETSMYKYPDSISYMPSYHAHQKQVNFVQPHPSMPVTSSSFFGAAS  
 QYWTSPTAGIYPNPSVPRHPNTHVPSHLGSYY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



<b>ACCN:</b>	NM_001017381
<b>ORF Size:</b>	1356 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001017381.1</a> , <a href="#">NP_001017381.1</a>
<b>RefSeq Size:</b>	1871 bp
<b>RefSeq ORF:</b>	1359 bp
<b>Locus ID:</b>	315532
<b>Cytogenetics:</b>	8q21
<b>MW:</b>	51 kDa