

Product datasheet for **RR203430L3V**

Slc17a6 (NM_053427) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc17a6 (NM_053427) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Slc17a6
Synonyms:	Dnpi; Vglut2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_053427
ORF Size:	1746 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR203430).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053427.1 , NP_445879.1
RefSeq Size:	3982 bp
RefSeq ORF:	1749 bp
Locus ID:	84487
UniProt ID:	Q9J112
Cytogenetics:	1q22
Gene Summary:	Na ⁺ -dependent inorganic phosphate cotransporter; involved in Ca ²⁺ -dependent glutamate release from astrocytes [RGD, Feb 2006]



[View online »](#)