

## Product datasheet for **RR203459**

### C1rl (NM\_001002804) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** C1rl (NM\_001002804) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** C1rl  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR203459 representing NM\_001002804  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCCGATCGCC

ATGTGCTGGTTGCTCCTTTGGGAACTCTCCATACTGCCACACAGGCCTCTGTGCTTTAGCCAGC  
 AGTTTCTCAGCAACTGACATCCCCGGATACCCAGAGCCATATCAAAGGCCAGGAGCCATGTCTGA  
 CATCGAGGCCCGGAGGGTTTGCCGTGAGACTGATCTCCAGGACTTTGATCTGGAGCCGTCGCCAGGC  
 TGTGAAGGCGACTCTGTCAACATTTCAACCAGAGGAACAGATGCAACCCGGCTGTGTGGCAGCAAGGTT  
 CCTCTCTGGGAGCCCTCCAATCAAATGGAGTTTGATCCTCCGGGAGGAGCCTTCGGTCTGACCTCCG  
 GGCACACTCCTCAAGAACAAGTACACACCTCCACAAAGGCTTCTGGCTCTCTACCAAGCCGCTGTG  
 AGCCAACCCAATGGGACGCTGAGGCTTCCACACACTGGAGCCAACCCCTCCGAGATTCAGAACCCT  
 GCCCAGGTCCTATTATAAGGAAGAACAGACAGGGACACTTAGCTGCCATCCTCTCGGAAGTGAAGGA  
 CAGACAGCGTGGAGAAGAGGTTCTGAGTGTGCGGCTGTGGAAGACCAGTTGTCCCATTGCCGAG  
 AACCCAAACACTTTTGGTTCTTCCAGAGCCAAGCCTGGCAACTCCCCTGGCAAGCCTTCACCAGTATCT  
 ATGGCCGTGGGGCGGAGCCCTGCTGGGTGACAGATGGATCCTGACGGCTGCCCATACCATTTCCCAA  
 GGATAGCATCTATCTCAGAAAGAACAAGACTGTGAATGTGTTTCTGGGTACACAGACGTAGACGAGCTG  
 CTGAAACTGGGGAACCATCCTGTCCGAAGGTTGGTGGTTACCCAGACTACCGCCAGAGGAGTCCCACA  
 ACTTTGATGGGGATATCGCCCTCCTGGAGCTTGAGCACAGGGTCCCTCTGGGGCCAGTCTCCTCCCGGT  
 CTGTCTGCCGACAACGAGACCTTACCACAGCGGCTGTGGGTTACATCAGTGGGTTCCGGTGTGGAG  
 ATGGGCTGGTTGACGACAAAGTTGAAATACTCAAAGCTGCCTGTAGCTCCCCGGGAGGCCTGTGAAGCCT  
 GGCTCCGCCAGAGGCAGCGAAGTGGAGTGTCTGACAATATGTTCTGTGTCGGGAAGAGATGCAGGT  
 CAACAGTGTCTGCCAGGGAGACAGTGGCAGCGTCTACGTGGTGTGGGACGACCGTCTCTCGCTGGGTG  
 GCCACGGGCATTGTATCCTGGGGCGTGGGTGTGGCAAGGGTACGGCTTCTACACCAAAGTGTCTCAGCT  
 ACGTGGATTGGATTAAGGGGTGATTGAATGTAAGACCGATGCCCGAGGCC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RR203459 representing NM\_001002804  
 Red=Cloning site Green=Tags(s)


MCWLLWLGILHTCPTQASVLLAQQFPQQLTSPGYPEPYIKGQESHADIEAPEGFVRLIFQDFDLEPSPG  
 CEGDSVTISTRGTDAIRLCGQQGSSLGSPPNQMEFVSSGRSLRLTFRAHSSKNKVTHLHKGFALYQAAV  
 SQPNGDAEAFITPGANPPEIQNHCPGYPYKKEEQGTLSCPSSRKWKDRQRGEEVPECVPCGRPVVPIAE  
 NPNTFGSSRAKPGNFPWQAFSTIYGRGGGALLGDRWILTAHTIFPKDSIYLKRNKTVNVFLGHTDVEDL  
 LKLGNHVPVRRVVHPDYRQEESHNFDDIALLELEHRVPLGPSLLPVCLPDNETLYHSGLWGYISGFGVE  
 MGWLTTKLKYSKLPVAPREACEAWLRQRQRTEVFSNMFCVGEEMQVNSVCQGDGSGSVYVWDDRALRW  
 ATGIVSWGVCCKGKGYFYTKVLSYVDWIKGVIECKDRCPEA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



<b>ACCN:</b>	NM_001002804
<b>ORF Size:</b>	1383 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001002804.1</a> , <a href="#">NP_001002804.1</a>
<b>RefSeq Size:</b>	1386 bp
<b>RefSeq ORF:</b>	1386 bp
<b>Locus ID:</b>	408246
<b>UniProt ID:</b>	<a href="#">Q6IE64</a>
<b>Cytogenetics:</b>	4q42
<b>MW:</b>	51.1 kDa
<b>Gene Summary:</b>	human homolog catalyzes the cleavage of the haptoglobin precursor protein [RGD, Feb 2006]