

Product datasheet for **RR203472**

Paqr4 (NM_001017377) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
 Product Name: Paqr4 (NM_001017377) Rat Tagged ORF Clone
 Tag: Myc-DDK
 Symbol: Paqr4
 Synonyms: MGC109357
 Vector: pCMV6-Entry (PS100001)
 E. coli Selection: Kanamycin (25 ug/mL)
 Cell Selection: Neomycin
 ORF Nucleotide Sequence: >RR203472 representing NM_001017377
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGCGTTCTGACCGGGCCGCGTCTTCTGGACTGGGCTAGCTCGCCGCCACACCTGCAGTTCAATAAGT
 TCGTATTGACCGGCTACCGCCGGCCAGCAGCGGCTCGGGCTGCCTGCGCAGTCTTTTCTACCTACACAA
 CGAGCTGGGCAACATCTACACACAGGGCTAGCCCTGCTGGGCTTCTGGTGTAGTGCCAAATGACCATG
 CCCTGGAGTCAGCTGGGCAAGGATGGCTGGCTAGGAGGCACACATTGTGTGGCCTGCCTGGCGCCCTG
 CAGCCTCTGTGCTCTATCACCTTTCATGTGCCACCAGGAGGCAGTCCCTGTGTACACCCGGCTCCTTGC
 CCTGGATATGTGTGGAGTCTGTCTTGTCAACACCTTGGGGCCCTACCCATCATCCACTGCACTCTGGCT
 TGCAGACCCTGGCTTCGCCAGCTGCCCTGATGGGTTACTGCACTGTGAGGTGTAGCCGGCTGGAGAG
 CCCTCACTGCCCCCTCCACCAGTGCCCGGCTTCGAGCCTTTGGTTGGCAAGCTGGGGCCCGCTGCTGGT
 GTTTGGGGCCCGTGGAGTGGGGCTGGGCTCAGGGGCTCCAGGCTCCCTGCCCTGCTACCTGCGCATGGAC
 GCACTGGCTCTGCTGGGAGGCTGGTGAATGTGGCAGCCTGCCAGAGCGCTGGGGCCCTGGCCGCTTCG
 ACTACTGGGGCAACTCCACCAAATCATGCACTTGCTGAGCGTGGGCTCCATCCTCCAGTCCATGCTGG
 GGTTGTGCTGACCTGCTCTGGGCTGCGGCCATGCCTGTCCCCAGAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001017377.1 , NP_001017377.1
RefSeq Size:	2033 bp
RefSeq ORF:	822 bp
Locus ID:	302967
Cytogenetics:	10q12
MW:	29.2 kDa