

Product datasheet for **RR203484L3V**

Fabp12 (NM_001134614) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fabp12 (NM_001134614) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Fabp12
Synonyms:	RGD1565000
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001134614
ORF Size:	396 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR203484).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001134614.1 , NP_001128086.1
RefSeq Size:	747 bp
RefSeq ORF:	399 bp
Locus ID:	499570
UniProt ID:	B7SUM8
Cytogenetics:	2q23
Gene Summary:	May play a role in lipid transport.[UniProtKB/Swiss-Prot Function]



[View online »](#)