

## Product datasheet for **RR203508L4V**

### Atp5g3 (Atp5mc3) (NM\_053756) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Atp5g3 (Atp5mc3) (NM_053756) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Atp5mc3
Synonyms:	Atp5g3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_053756
ORF Size:	426 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR203508).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_053756.1</a> , <a href="#">NP_446208.1</a>
RefSeq Size:	723 bp
RefSeq ORF:	429 bp
Locus ID:	114630
UniProt ID:	<a href="#">Q71S46</a>
Cytogenetics:	3q23
Gene Summary:	pancreatic expression is induced after ethanol consumption [RGD, Feb 2006]



[View online »](#)