

## Product datasheet for **RR203566L4V**

### Cyp26a1 (NM\_130408) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Cyp26a1 (NM_130408) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cyp26a1
Synonyms:	Cyp26
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_130408
ORF Size:	1491 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR203566).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_130408.2</a> , <a href="#">NP_569092.2</a>
RefSeq Size:	1754 bp
RefSeq ORF:	1494 bp
Locus ID:	154985
Cytogenetics:	1q53
Gene Summary:	member of the cytochrome P450 superfamily of monooxygenase enzymes; this gene encodes an endoplasmic reticulum-specific protein which acts on retinoids, including all-trans-retinoic acid (Ra) [RGD, Feb 2006]



[View online »](#)