

## Product datasheet for **RR203611L3V**

### Taok1 (NM\_173327) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Taok1 (NM_173327) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Taok1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_173327
ORF Size:	3003 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR203611).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_173327.1</a> , <a href="#">NP_775449.1</a>
RefSeq Size:	3312 bp
RefSeq ORF:	3006 bp
Locus ID:	286993
UniProt ID:	<a href="#">O88664</a>
Cytogenetics:	10q24
Gene Summary:	protein kinase; involved in regulation of the p38-containing stress-responsive MAP kinase pathway [RGD, Feb 2006]



[View online »](#)