

Product datasheet for **RR203622L4V**

Wdr5 (NM_001039034) Rat Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Wdr5 (NM_001039034) Rat Tagged ORF Clone Lentiviral Particle |
| Symbol: | Wdr5 |
| Synonyms: | MGC114572 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_001039034 |
| ORF Size: | 1002 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RR203622). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001039034.1 , NP_001034123.1 |
| RefSeq Size: | 3140 bp |
| RefSeq ORF: | 1005 bp |
| Locus ID: | 362093 |
| UniProt ID: | Q498M4 |
| Cytogenetics: | 3p12 |



[View online »](#)

Gene Summary:

Contributes to histone modification. May position the N-terminus of histone H3 for efficient trimethylation at 'Lys-4'. As part of the MLL1/MLL complex it is involved in methylation and dimethylation at 'Lys-4' of histone H3. H3 'Lys-4' methylation represents a specific tag for epigenetic transcriptional activation. As part of the NSL complex it may be involved in acetylation of nucleosomal histone H4 on several lysine residues. May regulate osteoblasts differentiation (By similarity).[UniProtKB/Swiss-Prot Function]