

Product datasheet for **RR203632**

Scfd2 (NM_001017499) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Scfd2 (NM_001017499) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Scfd2
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)


ORF Nucleotide
Sequence:

>RR203632 representing NM_001017499
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAGTGCCGTGGGCGTACTGTCTTTTGCCAGCAAGGCTGGGAGCAGGTGCTGGCCAAAGTGAAGTGGT
CCGTGGTCTACCTGGACGCTGCCTGCGCGGAGAGCCTGCACTGGAGCTGCGGTTCAAGCAGGCTTCTGGA
GGCGGTGAAGGGTCTGCGTGCAACCTGAGGGAATTTGAACCCCATGCAATAGGCGGTGGGGCCAAGCAG
CCCCGGGCGGTGTTCTGTTTGGAGCAGCCTGCTGAAAGGCCGCATAGTAGACACTCTACAAGACATCATT
GCCGCAGTCACTTTCAGCACTGTGTGGTGGTACAGCCGTGAGCCACCCAGTCCACCTCACTGCTAACCA
TGTACCAGCAGCTGCAGCGGCCGAGCTAGAGGGGCAGCAGCCGGTGTTCGAACAGCTGGAGGAGAAGCTC
TGTGAGTGGATGGCAACATGAACTACACAGCTGAGGTGCTCCACGTCCCCTTGTGCTCGCCCCTGCTG
CTCCCCACCTTGCTTTACTCCAGCCTTTGCTACCCTTTTCCACTGCTACCTCAGGATGTGCATGTTCT
TAATAATGCCAGACCCGACAAGAGGAGACTGAACAGCTTGGGTGAGGTGGATGCCACTGCCCTAACCCCA
GAGCTGCTGCTGTACATTAGATGCCTCGTGTGAGGCTCAGTTCTCTGTGTGAACATCTAGGGGTACGAG
AGGAGTGTTCGCTGTGGGTTCCCTCAGTCGGGTATCGCTACAGATCTGGCAAATTATGCCCTGCCAA
GAACAGGAGAAAGACTGCTACAGGCAGGGCTTCCGTGGTCTTTGTGGATAGAACTCTGGATCTCACGGGA
GCGGTTGGACACCATGGAGACAATTAGTAGAGAAGATCATTTCAGCACTTCCACAGCTTCCGGGCCACA
CCAACGACGTGATGGTCAACATGGCAGAGCTCACGGCTGTCCAGGCTGTGGAGGAAAACAGACCGTGAT
TGCACCAGGCTGCCTTGCCCATCCAATGATGCCTCAGCCAAGGCCCTCTGGGAAGCCTTACTGAACACC
AAGCACAAAGAGGCGGTGATGGAAGTCCGAGACATCTAGTGAAGCAGCCAGCAGAGAGAACCTGCCCA
TCAAGATGAGTATGGGGAGAGTCACTCCAGGGCAGCTCATGTCCATATTCAGCTTTCAAGAACAACCT
CAGAGCTTTACAAATCATTGTGGCTTCTACAGCTAGGGATGGCCACAGTTCAAACCTGAAGCACCCA
CAGACTGCCAAGTGGGATAAATTCTGGCATTGAAAGGCTTCTTCTACAGAGTCTTGGGATTCAACCA
TGGCTGGTGTGCTAAATCAGCTGCTACCCTTGATTAAGTCTCAAGCCAGAGAACCAGGGACGACTACAG
CCCCGAGGAGCTGCTAATCCTCCTCATATACATCTACTCGGTCACTGGGAATGCACGGGGACACAGAC
CTGGGAGACGTGGAAGAGAAGGTGAGGAGAGCACTGGCCAGTGTCTTCTGAGGAGTCTGAGCTGTCCC
CTTTGCTGCAGAAGATCACAGGCTGTGACTCTGCAATCGACCTGACATTGCCAAATCCCAAACCTGCTGT
GAATGATGTCTTACGGCTCTCCGGGAGATCGCTGGAGCTAGGAATCTCCTGAGACAGTTCAAATCGGTG
TATGTTCTGGAATAGTACCCATCAGGCAGCCTATAAGCCACTGTTGAAGCAGGTTGTGGAAGAAATAT
TTAATCCTGAGAAGCCAGATCCCATTGACATTGAACATATGTCCTCAGGCCTCACTGATCTCCTTAAAC
TGGATTACAGCATGTTTATGAAGGTGAGCCGGCCCCATCCCAGTGACCATCCCCTCTTGATCCTTTGTG
GTGGGCGGCGTCACAGTTGCTGAAGCCAAAATGGTCAAAGACCTTGTAGCATCTCTGAAGCCAGGAACCC
AGGTGATGGTGTGTCCACAAGACTCCTGAAGCCACTTAACATTCCTGAGCTGCTCTTGAACCTGACCC
GCTACACCCAGACCTTGGCTTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Plasmid Map:


ACCN: NM_001017499

ORF Size: 2052 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001017499.1](#), [NP_001017499.1](#)

RefSeq Size: 2360 bp
RefSeq ORF: 2055 bp
Locus ID: 498353
Cytogenetics: 14p11
MW: 74.8 kDa