

## Product datasheet for **RR203642L3V**

### Uba3 (NM\_057205) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Uba3 (NM_057205) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Uba3
Synonyms:	Ube1c
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_057205
ORF Size:	1386 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR203642).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_057205.2</a> , <a href="#">NP_476553.1</a>
RefSeq Size:	2189 bp
RefSeq ORF:	1389 bp
Locus ID:	117553
UniProt ID:	<a href="#">Q99MI7</a>
Cytogenetics:	4q34
Gene Summary:	catalytic subunit of the activating enzyme of the ubiquitin-like NEDD8 conjugation pathway (known as neddylation); inhibits transcription induced by steroid hormone receptors [RGD, Feb 2006]



[View online »](#)