

Product datasheet for **RR203658**

Actr10 (NM_001009602) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Actr10 (NM_001009602) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Actr10
Synonyms: MGC95127
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR203658 representing NM_001009602
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCCGCTCTATGAGGGCCTCGGGAGCGGGCGGAGAGAAGACGGCGGTCTGATCGACCTAGGAGAAGCCT
TCACCAAGTGTGGATTCGCAGGAGAACTGGTCCACGATGCATAATTCCTAGTGTGATAAAAAGAGCTGG
CATGTCTAAGCCGATCAAAGTTGTCCAGTATAATCAATACAGAAGAACTATATTCCTACCTAAAGGAA
TTCATCCACATACTGTATTTTCAGGCATCTGTTGGTGAATCCAGAGACCGCCGGTTGTGGTTATTGAGT
CAGTGTTATGTCCTTCCCACTTCAGAGAGACCCTTACTCGTGTTCTTTCAAATATTTTGGAGTCCCATC
TGTCTACTTGCTCCAAGTCATCTGATGGCACTCCTGACGCTCGGGATTAATTCGCCATGGTCTGGAT
TGTGGATATAGGGAAAAGCCTGGTGTTCCTATCTATGAAGGCATCCAGTACTGAATTGCTGGGGAGCTC
TGCCCTTAGGAGGAAAAGCTCTTCAAGGAGTTGAAACTCAGCTGTTGGAAACAATGCACTGTTGACAC
TGGTGTGCTAAAGGACAGAGCCTTCCCTCGGTGATGGTTTCAGTCCAGAAAGCGTCTAGAAGATATT
AAAGTCCGGACCTGCTTTGTAAGTGACCTGCAGCGTGGACTACAAATCCAAGCAGCAAAGTTAATTG
ATGGAAATAATGAGCGTCCCTCTCCACCTCCAAATGTCGACTATCCATTGGATGGAGAGAAAATTTACA
TGTGCTTGGATCAATCAGAGATTCAGTCGTGAAATTCCTTTTGAACAAGATAATGAGGAGAAGTCAGTT
GCCACTTTAATATTGGATTCCCTTTTACAGTGCCAGTAGACACCAGGAAGCAACTGGCAGAGAATCTGG
TAATCATAGGAGGCACGTCTATGTTGCCAGGATTTCTCCACAGATTGCTTGCAGAAATACGGTATTTGGT
AGAAAAGCCAAAATACAAAAGACACTTGGCACCAAGAAGTTCGAATTCATACTCCACCTGCAAAGGCT
AATTGTGTGGCTGGTTGGGAGGCGCTGTTTTGGAGCATTGCAAGACATACTGGGAGCCGCTCCATTT
CAAAGGAGTACTATAATCAGACGGGCCGACACCTGACTGGTGTCTCTGAACAACCCGCTCTGGAGAT
GATGTTTGATGTTGGGAAGACGCAGCCACCGCTGATGAAAAGAGCATTTTCCACCGAGAAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >RR203658 representing NM_001009602
 Red=Cloning site Green=Tags(s)


MPLYEGLGSGGEKTAVVIDLGEAFTKCGFAGETGPRCIIPSVIKRAGMSKPIKVVQYNINTEELYSYLKE
 FIHILYFRHLLVNPDRRVVIVIESVLCPSHFRETLTRVLFKYFEVPSVLLAPSHLMALLTLGINSAMVLD
 CGYRESLVLPYIEGIPVLNCWGALPLGGKALHKELETQLLEQCTVDTGAAKGQSLPSVMGSPVESVLEDI
 KVRTCFVSDLQRGLQIQAAKFNIDGNNERPSPPPNVVDYPLDGEKILHVLGSIIRDSVVEILFEQDNEEKSV
 ATLILDSELLQCPVDTRKQLAENLVIIGTSMPLPGFLHRLLAEIRYLVEKPKYKKTGTFKNFRIHTPPAKA
 NCVAWLGGAVFGALQDILGSRSSISKEYYNQTRTPDWCSLNNPPEMMFDVVGKTQPPLMKRAFSTEK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001009602

ORF Size:	1251 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001009602.1 , NP_001009602.1
RefSeq Size:	1457 bp
RefSeq ORF:	1254 bp
Locus ID:	299121
Cytogenetics:	6q24
MW:	46.2 kDa