

Product datasheet for **RR203671**

Pgrmc1 (NM_021766) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
 Product Name: Pgrmc1 (NM_021766) Rat Tagged ORF Clone
 Tag: Myc-DDK
 Symbol: Pgrmc1
 Synonyms: 25-Dx; 25Dx; MPR; VEMA
 Vector: pCMV6-Entry (PS100001)
 E. coli Selection: Kanamycin (25 ug/mL)
 Cell Selection: Neomycin
 ORF Nucleotide Sequence: >RR203671 representing NM_021766
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCTGCCGAGGATGTGGTGGCGACTGGCGCCGACCCAGCGAGCTGGAGGGCGGGGGCTGCTTCAAG
 AGATTTTCACGTCGCCTCTCAACCTGCTGCTCCTTGGCCTCTGCATCTTCTGCTCTACAAGATCGTTTCG
 CGGGGACCAGCCCGGTGCCAGTGGGGACAACGACGACGACGAGCCGCCCGCTGCCTCGCCTCAAGCCG
 CGTGACTTCACCCCTGCCGAATAAGGCGATACGATGGAGTCCAGGACCCGCGCATTCTTATGGCCATCA
 ACGGCAAGGTGTTTCGACGTGACCAAGGCCGCAAGTTCTATGGCCGGAGGGACCATACGGGGTCTTTGC
 TGGAAAGAGATGCATCCAGGGGCCTTGCCACATTTTGCCTGGACAAAGAAGCACTGAAGGATGAGTATGAT
 GACCTTTCTGACCTCACTCCTGCCAGCAGGAGACCCTGAATGACTGGGACTCTCAGTTCAGTTCACCTT
 CAAGTACCATCACGTGGGAAAAGTCTTGAAGGAGCGGAGGAGCCGATTGTGTACTCGGATGATGAAGA
 ACAAAGATGAGGCTGCTCGGAAGAGTGACTGAAGCAGTCAGTGGAGCATATCTATTTTTGTATTTTGCA
 AAATCATTGTAACTCCAGTCTGTCTTTACAACATGG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

| | |
|-------------------------------|---|
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_021766.1 , NP_068534.1 |
| RefSeq Size: | 1885 bp |
| RefSeq ORF: | 672 bp |
| Locus ID: | 291948 |
| UniProt ID: | P70580 |
| Cytogenetics: | Xq35 |
| MW: | 24.7 kDa |
| Gene Summary: | may be involved in regulation of axon guidance in the central nervous system [RGD, Feb 2006] |