

## Product datasheet for **RR203693L3V**

### Grk3 (NM\_012897) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Grk3 (NM_012897) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Grk3
Synonyms:	Adrbk2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_012897
ORF Size:	2064 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR203693).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_012897.2</a> , <a href="#">NP_037029.1</a>
RefSeq Size:	2373 bp
RefSeq ORF:	2067 bp
Locus ID:	25372
UniProt ID:	<a href="#">P26819</a>
Cytogenetics:	12q16
Gene Summary:	catalyzes phosphorylation of the beta-adrenergic receptor; may mediate desensitization of synaptic receptors; may regulate bone formation [RGD, Feb 2006]



[View online »](#)