

## Product datasheet for **RR203779**

### **Klkb1 (NM\_012725) Rat Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Klkb1 (NM_012725) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Klkb1
Synonyms:	Kalp1; KALP15; Klk3; PK
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RR203779 representing NM\_012725  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGATTTTATTCAAACAAGTGGGTATTTTGTTCCTTGTTCGCTACAGTTTCCTGTGGGTGTCTGTCCAC  
 AACTGTATGCAAATACCTTCTTCAGAGGTGGGATCTGGCTGCCATCTACACCCCGATGCCCAGCACTG  
 TCAGAAGATGTGCACGTTTCACCCAGGTGCCTGCTCTTCAGCTTCTTGCCGTGAGTCCAACCAAGGAG  
 ACAGATAAAAGGTTGGGTGCTTCATGAAAGAGAGCATTACAGGGACTTTGCCAAGAATACACCGGACAG  
 GGGCCATTTCTGGTCATTCTTTAAACAGTGTGGCCATCAATTAAGTGTGGCCACCAAGACATATACGA  
 AGGACTGGATATGAGAGGGTCCAACCTTAAATATATCTAAGACCGACAGTATTGAAGAATGCCAGAAACTG  
 TGCACAAATAATATTCAGTCCAAATTTTACATATGCTACAAAAGCATTTCACAGACCAGAGTACAGGA  
 AGAGTTGCCTGTGAAGCGCAGTTCAAGTGAACGCCACCAGTATAAAGCCAGTGGACAACCTGGTGTCT  
 TGGATTCTCACTGAAGTCTGTGCTCTCTCAGAGATCGGTTGCCCATGGATATTTCCAGCACTTTGCC  
 TTTGCAGACCTGAATGTAAGCCATGTGCTCACCCCGATGCCTTCGTGTGTCGCACCGTTTGCACCTTCC  
 ATCCCAACTGCCTCTTCTTACATTCTACACGAATGAGTGGGAGACGGAATCACAGAGGAATGTTTGTTT  
 TCTTAAGACATCTAAAAGTGAAGACCAAGTCCCCCTATTATTCAAGAAAATGCTGTATCTGGATACAGT  
 CTCTTACCTGCAGAAAAGCTCGCCCTGAACCTGCCATTTCAAGATTTACTCTGGAGTTGCCTTCGAAG  
 GGAAGAAGTGAACGCGACCTTCGTGCAGGGAGCAGATGCGTGCCAAGAGACTTGTACAAAGACCATCCG  
 CTGTCAAGTTTTTACTTACTCATTGCTTCCCAAGACTGCAAGGCAGAGGGGTGTAATGTTCTTAAAGG  
 TTATCCACGGATGGCTCTCCAAGTAGGATCACCTATGAGGCACAGGGGAGCTCTGGTTATTCTCTGAGAC  
 TGTGTAAGTTGTGGAGAGCTCTGACTGTACGACAAAAATAAATGCACGTATTGTGGGAGGAACAACTC  
 TTCTTTAGGAGAGTGCCATGGCAGGTCAGCCTGCAAGTAAAGTTGGTTTCTCAGAATCATATGTGTGGA  
 GGGTCCATCATTGGACGCCAATGGATACTGACGGCTGCCATTGCTTTGATGGGATTCCTATCCAGACG  
 TGTGGCGTATATATGGCGGGATTCTTAATCTGTGAGAGATTACAAACAAAACGCCTTTCTCAAGTATAAA  
 GGAGCTTATTATTCATCAGAAATACAAAATGTCAGAAGGCAGTTACGATATTGCCTTAATAAAGCTTCAG  
 ACACCGTTGAATTATACTGAATTCAAAAACCAATATGCCTGCCTTCCAAAGCTGACACAAATACAATTT  
 ATACCAACTGCTGGGTGACTGGATGGGGCTACACAAAGGAACGAGGTGAGACCCAAAATATTCTACAAA  
 GGCAACTATTCCCTTGGTACCAATGAAGAATGCCAGAAAAAATATAGAGATTATGTTATAACCAAGCAG  
 ATGATCTGTGCTGGCTACAAAGAAGGTGGAATAGATGCTTGTAAAGGAGATTCCGGTGGCCCTTAGTTT  
 GCAAACATAGTGAAGGTGGCAGTTGGTGGGTATCACCAGCTGGGGCGAAGGCTGTGCCCGCAAGGAGCA  
 ACCAGGAGTCTACACAAAGTTGCTGAGTACATTGACTGGATATTGGAGAAGATACAGAGCAGCAAGGAA  
 AGAGCTCTGGAGACATCTCCAGCA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR203779 representing NM\_012725  
Red=Cloning site Green=Tags(s)


```
MILFKQVGYFVSLFATVSCGCLS QLYANTFFRGGDLAAIYTPDAQHCQKMCTFHPRCLLFSFLAVSPTKE  
TDKRFGCFMKESITGTLPRIHRTGAI SGHSLKQCGHQLSACHQDIYEGLDMRGSNFNISKTDSEIECQKL  
CTNNIHCQFFTYATKAFHRPEYRK SCLLKRSSSGTPTSIKPVDNLVSGFSLKSCALSEIGCPMDIFQHFA  
FADLNVSHVVPDAFVCRVTCTFHPNCLFFTFYTNNEWETESQRNVCFLKTSKSGRSPPIIQENAVSGYS  
LFTCRKARPEPCHFKIYSGVAFEGEELNATFVQGADACQETCTKTIRCQFFTYSLLPQDCKAEGCKCSLR  
LSTDGSPTRITYEAQGSSGYSRLCKVVESSDCTTKINARIVGGTNSLGEWPQVSLQVKLVSNHMC G  
GSIIGRQWILTAAHCFDGIPYPDVWRIYGGILNLSEITNKTPFSSIKELIIHQYKMSEGSYDIALIKLQ  
TPLNYTEFQKPICLPSKADTNTIYTNWVTGWGYTKERGETQNILQKAT IPLVPNEECQKKYRDYVITKQ  
MICAGYKEGGIDACKGDSGGPLVCKHSGRWQLVGITSWGEGCARKEQPGVYTKVAEYIDWILEKIQSSKE  
RALETSPA
```

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/ja2016\\_c03.zip](https://cdn.origene.com/chromatograms/ja2016_c03.zip)

**Restriction Sites:** Sgfl-MluI

**Cloning Scheme:**



**ACCN:** NM\_012725

**ORF Size:** 1914 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_012725.2](#), [NP\\_036857.2](#)

**RefSeq Size:** 2513 bp

**RefSeq ORF:** 1917 bp


**Locus ID:** 25048

**Cytogenetics:** 16q11

**MW:** 71.7 kDa

**Gene Summary:** serine protease that plays a role in blood coagulation [RGD, Feb 2006]

Product images:



Circular map for RR203779