

## Product datasheet for **RR203880L3V**

### Pip4k2b (NM\_053550) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Pip4k2b (NM_053550) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Pip4k2b
Synonyms:	Pip5k2b
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_053550
ORF Size:	1248 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR203880).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_053550.1</a> , <a href="#">NP_446002.1</a>
RefSeq Size:	1987 bp
RefSeq ORF:	1251 bp
Locus ID:	89812
UniProt ID:	<a href="#">O88377</a>
Cytogenetics:	10q31
Gene Summary:	has similarity to phosphatidylinositol-5-phosphate 4-kinase IIalpha and IIgamma isoforms [RGD, Feb 2006]



[View online »](#)