

Product datasheet for **RR203886**

Taf9b (NM_133615) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Taf9b (NM_133615) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Taf9b
Synonyms: Dn7; Taf9l; TAFII31
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR203886 representing NM_133615
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGCCTCTCAAGAACGCTCCGAGGGATGCCTTGGTGATGGCACAGATCCTGAAGGATATGGAATCA
CAGATTATGAGCCAAGGGTTATAAATCAAATGTTGGAATTTGCCTCCGATATGTGACTACAATTCTGGA
TGATGCAAAAATTTATTCAAGTCAGCTAAGAAACCTACTGTGGATGCAGATGATGTGAGACTGGCCATC
CAGTGTCCGGGCTGACCAGTCTTTACCTCTCCTCCTCCAAGAGATTTTTACTGGATATTGCAAGGCAGA
AGAATCAAACCCCTCTGCCGCTGATTAACCATATGCAGGACCTAGATTGCCACCTGATAGATACTGCTT
AACAGCTCCAAACTATAGGCTGAAGTCCTTAGTTAAAAAGGGACCTAACCAAGGAAGACTAGTTCCTAGA
TTAAGTGTCCGGTCCGGTCCAGTAGCAGACCTACTACTCCTACTGTAGCACCCCAAGCAGTATCTGTCC
CAAATAAGGTTGCAACTCCAGTGTGAGTGAAGCAAGCAAGATTTGCCGTGCAGATTCCAACCTCTCAGTC
CACGCTGCCAAACCAGCTCCTGCAGCGACTTCAGTCCAAATGTTCTGATTAATCCTTCAATGATTGGG
TCCACAAATATTCTCATTACCACCAACATGGTTTCCTCACAGAGCACAGCCACTGACTCAAACCCACTGA
AGAGGAAACATGATGATGATGATGATGACGACGACGATGATGATGACAATGACACTATG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR203886 representing NM_133615
 Red=Cloning site Green=Tags(s)

MAPLKNAPRDALVMAQILKDMGITDYEPVINVQMLEFAFRVYVTTILDDAKIYSSHAKKPTVDADDVRLAI
 QCRADQSFTSPPPRDFLLDIARQKNQTLPLIKPYAGPRLPPDRYCLTAPNYRLKSLVKKGPNQGRVLR
 LSVGAVSSRPTTPTVAPPQAVSPNKVATPVSVTSQRFVAVQIPTSQSTPAKPAPAATSVQNVLINPSMIG
 STNILITTNMVSSQSTATDSNPLKRKHDDDDDDDDDDDDNDMT


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_133615
ORF Size: 759 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_133615.1 , NP_598299.1
RefSeq Size:	2477 bp
RefSeq ORF:	762 bp
Locus ID:	171152
UniProt ID:	Q62880
Cytogenetics:	Xq22
MW:	27.6 kDa
Gene Summary:	clone which is induced in differentiated PC12 cells during apoptosis; transcriptional coactivator of Tp53 protein [RGD, Feb 2006]