

Product datasheet for **RR203897**

Cd59 (NM_012925) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cd59 (NM_012925) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Cd59
Synonyms: Cd59a; Cd59b; MAC-IP; MACIF; MACIP
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR203897 representing NM_012925
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAGAGCTCGGAGGGGATTCATCTTACTCCTGCTTCTGGCTGTCCTCTGTTCCACAGGTGTTAGCCTCA
GATGCTACAACCTGTTTAGACCCGGTTTCTTCATGCAAAACAAACAGCACTTGCTCTCCTAACCTGGATGC
TTGTCTTGTGTGTATCCGAAAGCAAGTCTATCAACAGTGTTGGAGATTTTCGGATTGTAATGCCAAG
TTCATTTTGGCCGACTAGAAATCGAAACGTACAATACAGATGCTGCCAGGCGGACTTGTGTAACAAA
GCTTCGAAGACAAGCAAACAATGGGGCAATCTCCCTATTGGGAAGACAGCGTTGCTGGTGACCTCGGT
TCTGGCGGCCATTTGAAGCCTTGTTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR203897 representing NM_012925
Red=Cloning site **Green**=Tags(s)

MRARRGFILLLLAVLCSTGVSLRCYNCLDPVSSCKTNSTCSPNLDACLVAVSGKQVYQQCWRFSDCNAK
FILSRLEIANVQYRCCQADLCNKSFEKPNNGAISLLGKTALLVTSVLAAILKPCF


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_012925

ORF Size: 378 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_012925.1 , NP_037057.1
RefSeq Size:	1523 bp
RefSeq ORF:	381 bp
Locus ID:	25407
UniProt ID:	P27274
Cytogenetics:	3q32
MW:	13.8 kDa
Gene Summary:	may play a role in regulation of complement; may bind complement proteins C8 and C9 [RGD, Feb 2006]