

Product datasheet for **RR203988L3V**

Tacr1 (NM_012667) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tacr1 (NM_012667) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Tacr1
Synonyms:	Tac1r
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_012667
ORF Size:	1221 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR203988).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012667.2 , NP_036799.1
RefSeq Size:	3408 bp
RefSeq ORF:	1224 bp
Locus ID:	24807
UniProt ID:	P14600
Cytogenetics:	4q34
Gene Summary:	mediates neuropeptide signaling; may play a role in pain sensitivity; may mediate baroreflex transmission for blood pressure regulation [RGD, Feb 2006]



[View online »](#)