

Product datasheet for **RR203989**

Rnase10 (NM_001012467) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rnase10 (NM_001012467) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Rnase10
Synonyms: Rah1; RGD1310415
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR203989 representing NM_001012467
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAAGGTGACACTGGTGCATCTGTTGTTTCATGATGTTACTGCTGTTGCTAGGACTGGGGTGGGGCTGG
GGCTGGGCCTTCACATGGCTGCTGCCATCCTGGAGAATCAGCCACTGGATGAGTTTTGGCCAGTGACTC
CCAGGACACAGCTGAGGCCACTGAAGAGGGACAGGGCACCAGGACCACAGAAGCCCTGGTACTTGACAAC
AAAGAAATGGCAGTACCTGTCTGGTCAGAAGATACTGTCCTCAGTGAAGATGAAGTGGGGGAAGCAGGA
TGCTGAGAGCCAAGACCCTCCTTCAGAGCAAGCAAGGCTACCTTAAGTTTGACTTGAACATCCGGGACTG
TAACGTGATGATGGCCACAAGATAAAGGAGCACAATCAGAGCTGCATAAACGATTATACATTCATCCAT
GAGGATCCAAGCACAGTCGGAGCCGTCTGTAACAGTCCCCTGGTTGACTGTGATCTCAAGGGAGGCAAA
GTCACAAAAGCCCCGGCCTTTTGACCTGACATTGTGCAAGTTGGCCAAACCAGGCCAAGTCACTCCCAA
CTGCCACTATCTGACTTACATAACTGAAAAGGTCATTATCATAACATGCAACAACACGAAACAACCTGGAG
ATTAAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

| | |
|-------------------------------|---|
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_001012467.1 , NP_001012485.1 |
| RefSeq Size: | 639 bp |
| RefSeq ORF: | 639 bp |
| Locus ID: | 305840 |
| UniProt ID: | Q5GAM0 |
| Cytogenetics: | 15p14 |
| MW: | 23.5 kDa |
| Gene Summary: | Secreted proximal epididymal protein required for post-testicular sperm maturation and male fertility. May be involved in sperm adhesion to the egg zona pellucida. Does not have ribonuclease activity (By similarity).[UniProtKB/Swiss-Prot Function] |