

## Product datasheet for RR204046L4

### Gosr2 (NM\_031685) Rat Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Gosr2 (NM_031685) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Gosr2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204046).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.



[View online »](#)



RefSeq ORF: 639 bp  
Locus ID: 64154  
UniProt ID: [O35165](#)  
Cytogenetics: 10q32.1  
Gene Summary: may play a role in vesicle transport from the cis/medial to the trans-Golgi/TGN [RGD, Feb 2006]