

Product datasheet for **RR204051L4V**

P2ry6 (NM_057124) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	P2ry6 (NM_057124) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	P2ry6
Synonyms:	P2y6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_057124
ORF Size:	984 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204051).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_057124.2 , NP_476465.1
RefSeq Size:	2530 bp
RefSeq ORF:	987 bp
Locus ID:	117264
UniProt ID:	Q63371
Cytogenetics:	1q32
Gene Summary:	purinergic nucleotide receptor that displays increased activation by UTP compared to ADP; activation induces phospholipase C activation and calcium mobilization [RGD, Feb 2006]



[View online »](#)