

## Product datasheet for **RR204077L3V**

### Xpnpep2 (NM\_057155) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Xpnpep2 (NM_057155) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Xpnpep2
Synonyms:	mAPP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_057155
ORF Size:	2022 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204077).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_057155.3</a> , <a href="#">NP_476496.1</a>
RefSeq Size:	2828 bp
RefSeq ORF:	2025 bp
Locus ID:	117522
UniProt ID:	<a href="#">Q99MA2</a>
Cytogenetics:	Xq36
Gene Summary:	human homolog is a candidate gene for premature ovarian failure [RGD, Feb 2006]



[View online »](#)