

Product datasheet for **RR204084**

Prss21 (NM_181477) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Prss21 (NM_181477) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Prss21
Synonyms: Esp-1; Resp-1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR204084 representing NM_181477
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAGCGCGGGGCAAGACGTTGGTGCCACTGCTGGTGGTTGTGGTGGCGGTGGAGGTGACCTTACAGT
CAACCTCTTCGCATGTGAAACCCGTGGATCCCGAGAAACCCGAAGTGCAGGAAGCCAACCTATTGTCAGG
ACCCTGCGGTACAGGACCATCCCTTCAGTATAGTGGGTGGCGAAGAGGCCGAGCTTGGCCGCTGGCCT
TGGCAAGGGAGCCTGCGTGTATGGGGCAACCCTTATGTGGCGCAACTTTGCTCAACCGCCGCTGGGTGC
TTACAGCTGCCCACTGCTTCCAAAAGGACAATGATCCTTTCGACTGGACAGTCCAGTTTGGCGAGCTGAC
CTCCAGGCCGTCTCTCTGGAACCTACAGGCTTATTCCAACCGTTACCAAAATAGAAGATATTTTCTGAGC
CCCAAGTACACGGAGCAGTTTCCCATGACATAGCTCTGCTGAAGCTGTCATCGCCTGTCACCTACAGTA
ACTTCATCCAGCCCATCTGCCTCCTGAACCTCCACATAAAGTTTGGCAACCGAAGTACTGCTGGGTGAC
CGGCTGGGGGGCTATTGGAGAGGATGAGAGTCTGCCACTCCCAACAACCTCCAGGAAGTGCAGGTAGCT
ATCATCAACAACCATGTGTAACCATTTGTTCAAAAAGCCAGACTTCCGAATAAACATCTGGGGAGACA
TGTTTTGCGCTGGCAGTCTGAAGCGGCAAGGATGCCTGCTTTGGTGACTCAGGAGGACCTTGGTCTG
CAACCAGGATACAGTGTGGTATCAGGTCGGAGTTGTGAGCTGGGAATAGGCTGTGGTCGGCCCAATCGC
CCTGGAGTCTATACCAACATCAGTCATCACTACAACCTGGATCCGTTGACTATGATCCGAAACGGGATGC
TCAGGCCGTGACCCAGCCCCCTACTGCTGTTTCTACTCTGTCCTGGCTTCTTCTCTCTGCTGAGGCC
TGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR204084 representing NM_181477
Red=Cloning site Green=Tags(s)

MSARGKTLVPLL VVVVAVEVTLQSTSSHVKPVDPEKPELQEANLLSGPCGHRTIPSRIVGGEEAELGRWP
 WQGLRVWGNHLCGATLLNRRWVLTAAHCFQKDNDFDWTVQFGLTSRPSLWNLQAYSNRYQIEDIFLS
 PKYTEQFPHDIALLLKSSPVTYSNFIQPICLLNSTYKFNRTDCWVTGWGAIGEDSLPLPNNLQEVQVA
 IINNTMCNHLFKKPDFRINIWGMVYCAGSPEGGKDACFGDSGGPLVCNQDTVWYQVGVVSWGIGCGRPNR
 PGVYTNISHHYNWIRLTMIRNGMLRDPAPLLLLFLTLAWASSLLRPA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_181477

ORF Size: 984 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_181477.2](#), [NP_852142.1](#)

RefSeq Size: 1101 bp

RefSeq ORF: 987 bp

Locus ID: 353251

Cytogenetics: 10q12

MW: 36.6 kDa

Gene Summary: serine protease; may be involved in spermatogenesis and fertilization [RGD, Feb 2006]