

Product datasheet for **RR204203L3V**

Crh (NM_031019) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Crh (NM_031019) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Crh
Synonyms:	CRF
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031019
ORF Size:	561 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204203).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031019.1 , NP_112281.1
RefSeq Size:	1197 bp
RefSeq ORF:	564 bp
Locus ID:	81648
UniProt ID:	P01143
Cytogenetics:	2q24
Gene Summary:	binds CRH receptors and mediates neuroendocrine stress response; may play a role in the consolidation of long term memory [RGD, Feb 2006]



[View online »](#)