

## Product datasheet for **RR204276L4V**

### Ywhaz (NM\_013011) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ywhaz (NM_013011) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ywhaz
Synonyms:	14-3-3z
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_013011
ORF Size:	735 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204276).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_013011.3</a> , <a href="#">NP_037143.2</a>
RefSeq Size:	1815 bp
RefSeq ORF:	738 bp
Locus ID:	25578
UniProt ID:	<a href="#">P63102</a>
Cytogenetics:	7q22
Gene Summary:	30kDa component of the mitochondrial import stimulation factor, a protein complex that facilitates the import of in vitro synthesized precursor proteins into mitochondria [RGD, Feb 2006]



[View online »](#)