

## Product datasheet for **RR204338L4V**

### Dad1 (NM\_138910) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Dad1 (NM_138910) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Dad1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_138910
ORF Size:	339 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204338).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_138910.3</a> , <a href="#">NP_620265.1</a>
RefSeq Size:	722 bp
RefSeq ORF:	342 bp
Locus ID:	192275
UniProt ID:	<a href="#">P61805</a>
Cytogenetics:	15p13
Gene Summary:	human homolog inhibits apoptosis [RGD, Feb 2006]



[View online »](#)