

## Product datasheet for **RR204347L4V**

### **lapp (NM\_012586) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	lapp (NM_012586) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	lapp
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_012586
ORF Size:	279 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204347).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_012586.2</a> , <a href="#">NP_036718.1</a>
RefSeq Size:	282 bp
RefSeq ORF:	282 bp
Locus ID:	24476
UniProt ID:	<a href="#">P12969</a>
Cytogenetics:	4q44
Gene Summary:	pancreatic islet hormone found in the pancreatic amyloid deposits [RGD, Feb 2006]



[View online »](#)