

Product datasheet for **RR204362L4V**

Coro6 (NM_139115) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Coro6 (NM_139115) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Coro6
Synonyms:	MGC112600
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_139115
ORF Size:	1416 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204362).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_139115.2 , NP_620815.2
RefSeq Size:	2444 bp
RefSeq ORF:	1419 bp
Locus ID:	245982
UniProt ID:	Q920J3
Cytogenetics:	10q24
Gene Summary:	contains a WD40 domain [RGD, Feb 2006]



[View online »](#)