

## Product datasheet for **RR204387L1V**

### **Etv4 (NM\_001108299) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Etv4 (NM_001108299) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Etv4
Synonyms:	Pea3
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001108299
ORF Size:	1458 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204387).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001108299.1</a>
RefSeq Size:	2329 bp
RefSeq ORF:	1461 bp
Locus ID:	360635
Cytogenetics:	10q32.1



[View online »](#)