

Product datasheet for **RR204396L3V**

Actr3 (NM_031068) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Actr3 (NM_031068) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Actr3
Synonyms:	Arp3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031068
ORF Size:	1254 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204396).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031068.1 , NP_112330.1
RefSeq Size:	2114 bp
RefSeq ORF:	1257 bp
Locus ID:	81732
UniProt ID:	Q4V7C7
Cytogenetics:	13q12
Gene Summary:	actin related protein that may participate in the organization of filament/membrane adhesions [RGD, Feb 2006]



[View online »](#)