

Product datasheet for **RR204489**

Papd5 (NM_001107416) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Papd5 (NM_001107416) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Papd5
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR204489 representing NM_001107416
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCTCCAGACCCGAGGAGGAGAAGATGCGGATGGAGGTGGTGAAGCAGGATCGAGAGTGTGATTAAGG
AGCTCTGGCCAGTGTGACTGGGAGAACCTGCCTCTTTGGACCCTGGAAGAAGCCCTTCGAAACACAA
AGTCGTGACGAGGATCCGTGAAAGTTCTAGACAAAGCAACGGTTCCTATTATTAATTAACAGATTCT
TTTACTGAAGTGAAGTTGATATCAGCTTAATGTCCAGAATGGCGTGAGAGCAGCAGACCTCATTAAAG
ATTTTACTAAGTTACATCCCAGGGAAGATGCTTGCATCCCAATACAACTATGGTGTCTCTTAATAGA
ATTTTTGAATTATGGACGGCACTTCAATTACTTAAAGACTGGCATCCGGATAAAGGATGGTGGTTCC
TATGTGGCCAAAGATGAAGTACAGAAAAACATGCTTGATGGCTACCGCCATCAATGCTTATATCGAAG
ATCCTTTACAGCCAGGTAATGATGTTGGGAGGAGTTCGTATGGGGCCATGCAGGTAAAGCAGGCCTTTGA
CTACGCCATGTGGTGTGAGCCACGCAGTGTGCGCGATTGCAAAGTGTATCCCAACAACGAGACAGAG
AGCATATTGGGTAGAATAATTAGAGTGACGGATGAAGTTACCACATACAGAGATTGGATATCAAAACAGT
GGGGCTTGCGAATAGGCCCGAGCCTTCATGCAATGGAAATGGTGTACCTTGATAGTAGAGACGCAGCA
GTTAGATAAGTGTAAATAAATCTGTCTGAAGAGAATGAAGCCCCGGGAAAATGTAGAAGTAAAGCCTCG
GAGTCTCTTAGTAAACACTCTTCAAACCTTTCATCAGGTCTGTGCTCCTCCTCTGCCACACAGTCCA
GCTCTAGTGACGTTGATTCTGACGCAACACCATGCAAAACCCCCAAACAAGTGTCTGCCGCCACCCAC
TGTGAGCCGGGTAGGCTCACAGGATGTCTCCTTGGAGGTCTCTCAGGCAGTTGGGAAAATGCAGAGCACC
CAAACCACCAACCGCCAACAGCACCACCAAGTACAGCATGGATCAGCAAGGCTCTCCGTTCTTCCA
GCAAAGGCTTCCAAGGTACAGCTCAAACAGCCATGGTGCCTTGATGACAAGCAAACAGCATCCAGGCAA
ATCCAATACTCAGTGTACCATGGCAAAAAGAGGAGACACAAGAGGGACGCGCCCTCTCAGAGCTTTGT
AGG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



ACCN:	NM_001107416
ORF Size:	1263 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001107416.1 , NP_001100886.1
RefSeq Size:	3979 bp
RefSeq ORF:	1266 bp
Locus ID:	307745
Cytogenetics:	19p11
MW:	46.8 kDa