

Product datasheet for **RR204510L4V**

Dynll2 (NM_080697) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dynll2 (NM_080697) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Dynll2
Synonyms:	Dlc2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_080697
ORF Size:	267 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204510).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_080697.2 , NP_542428.1
RefSeq Size:	2370 bp
RefSeq ORF:	270 bp
Locus ID:	140734
UniProt ID:	Q78P75
Cytogenetics:	10q26
Gene Summary:	gephyrin-binding protein that plays a role in postsynaptic targeting of neurotransmitter receptors or scaffolding proteins [RGD, Feb 2006]



[View online »](#)