

## Product datasheet for **RR204537L3V**

### Gpr4 (NM\_001025680) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Gpr4 (NM_001025680) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Gpr4
Synonyms:	MGC116369
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001025680
ORF Size:	1095 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204537).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001025680.1</a> , <a href="#">NP_001020851.1</a>
RefSeq Size:	1513 bp
RefSeq ORF:	1098 bp
Locus ID:	308408
UniProt ID:	<a href="#">Q4KLH9</a>
Cytogenetics:	1q21
Gene Summary:	Proton-sensing receptor coupled to several G-proteins, including G(s), G(13) and G(q)/G(11) proteins, leading to cAMP production.[UniProtKB/Swiss-Prot Function]



[View online »](#)