

## Product datasheet for **RR204571L3V**

### **Pde1c (NM\_031078) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Pde1c (NM_031078) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Pde1c
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031078
ORF Size:	2304 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204571).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_031078.2</a> , <a href="#">NP_112340.1</a>
RefSeq Size:	3186 bp
RefSeq ORF:	2307 bp
Locus ID:	81742
UniProt ID:	<a href="#">Q63421</a>
Cytogenetics:	4q24
Gene Summary:	displays Ca <sup>2+</sup> /calmodulin dependent cAMP and cGMP phosphodiesterase activity; may play a role in olfactory signal transduction [RGD, Feb 2006]



[View online »](#)