

## Product datasheet for **RR204572L4V**

### Slu7 (NM\_001100550) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Slu7 (NM_001100550) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Slu7
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001100550
ORF Size:	1758 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204572).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001100550.1</a> , <a href="#">NP_001094020.1</a>
RefSeq Size:	2256 bp
RefSeq ORF:	1761 bp
Locus ID:	303057
UniProt ID:	<a href="#">Q80ZG5</a>
Cytogenetics:	10q21
Gene Summary:	human homolog is a splicing factor required for selection of the correct AG dinucleotide at the 3 splice site in the second catalytic step of mRNA splicing [RGD, Feb 2006]



[View online »](#)