

Product datasheet for **RR204576**

Adipor2 (NM_001037979) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Adipor2 (NM_001037979) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Adipor2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR204576 representing NM_001037979
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGAACGAGCCAACCGAACCCGATTGGGGTGCACCAGGACTCCAGAGCCAGATATAAGGCTCAGAAAAG
GGCACCAACTTGATGATACAAGGGGAGGTAACAATGACAACCACCGAGATTTGGAGCCAGTTTAGA
GACACCTGTTTGTCTTCTACTATGAAAATAGCCCCGAGGAACCTCGAGTGTGATGACAACCTCCAA
GAAGATGAAGGGTTTATGGGCATGTCCCTCTCCTACAAGCCCATCATGCTATGGAACGAATGGAAGAGT
TTGTTTGAAGGTGTGGGAAGTTCGATGGCGAGTATCCCTCATGATGTGCTGCCAGATTGGCTTAAGGA
CAATGACTTCCTTCTCCATGGACACCGCCTCCTATGCCTTCTTCCGGGCTGTTTTAAGAGCATCTTT
AGAATACACACAGAGACGGGCAACATTTGGACACATCTCCTAGGTTGTGATTCTTCTGTGCTGGGGA
TCTTTTATATGTTTCGCCAAATATATCTTTGTGGCCCTCTGCAAGAGAAAGTGGTCTTCGGGTTGTT
CTTCTGGGAGCCATTCTGCCTTCTTTTCATGGCTCTTCCACCGGTACTGCCACTCAGAAGGG
GTCTCCGACTCTTCTAAATTGGATTACTCTGGTATTGCTCTTCTGATCATGGGAAGTTTGTTCCTT
GGCTTTATTATTCTTCTACTGTAACCCACAACCTTGCTTCATCTACCTGATTGTCATCTGTGTGCTGGG
CATTGCAGCCATTATCGTCTCTCAGTGGGACATGTTTGCCACCCCTCAGTATCGGGGGTTCAGAGCAGGA
GTGTTCTGTGGGCTTGGGTCTGAGTGGAATCATCCCTACCTTGCAATATGTCATCTCAGAAGGGTTCTGA
AGGCTGCCACCATAGGGCAGATAGGCTGGCTAATGCTCATGGCTAGCCTCTATATCACAGGAGCTGCCCT
CTACGGCGCCCGTATCCCTGAGCGCTTCTTCTGGCAAATGTGACATCTGGTTTCACTTTCATCAGCTC
TTCCACATCTTTGTGGTTGCTGGTGCTTTGTTCACTTCCATGGAGTCTCAAACCTGCAGGAATTCGTT
TCATGATTGGTGGGGCTGCACTGAAGAGGATGCGCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >RR204576 representing NM_001037979
 Red=Cloning site Green=Tags(s)

MNEPTEHRLGCTRTPEDIRLRKGHQLDDTRGGNNDNHHGDLEPSLETPVCSSYYENSPEELEDHDDNSQ
 EDEGFMGMSPLLQAHHAMERMEEFVCKVWEGRWVIPHDVLPDWLKDNDLHGHRRPMPSPFRACFKSIF
 RIHTETGNIWTHLLGCVFFLCLGIFYMFRPNISFVAPLQEKVVFLGAILCLSFSLFHTVYCHSEG
 VSRLFSLDYSGIALLIMGSFVPLYYSFYCNQPFCFIYLVIVICVLGIAAIIIVSQWDMFATPQYRGVRAG
 VFVGLGLSGIIPTLHYVISEGFLKAATIGQIGWMLMASLYITGAALYAARIPERFFPGKCDIWFHFHQL
 FHIFVVAGAFVHFHGVSNLQEFRFMIGGCTEEDAL


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001037979

ORF Size:	1158 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001037979.1 , NP_001033068.1
RefSeq Size:	1185 bp
RefSeq ORF:	1161 bp
Locus ID:	312670
Cytogenetics:	4q42
MW:	44 kDa