

## Product datasheet for **RR204595L3V**

### IL1a (NM\_017019) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	IL1a (NM_017019) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	IL1a
Synonyms:	IL-1 alpha; IL-1F1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_017019
ORF Size:	810 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204595).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_017019.1</a> , <a href="#">NP_058715.1</a>
RefSeq Size:	1992 bp
RefSeq ORF:	813 bp
Locus ID:	24493
UniProt ID:	<a href="#">P16598</a>
Cytogenetics:	3q36
Gene Summary:	proinflammatory cytokine [RGD, Feb 2006]



[View online »](#)