

Product datasheet for **RR204597**

Ctgf (NM_022266) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ctgf (NM_022266) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Ctgf
Synonyms: CTGRP
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR204597 representing NM_022266
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGCTCGCCTCCGTCGCGGGTCCCCTAGCCTCGCCTGGTGCTCCTCCTCTGCACCCGGCCTGCCACCG
GCCAGGACTGCAGCGCAGTGTCACTGCGCAGCTGAAGCGCGCCGCGCTGCCCGCCGGCGTGAGCCT
GGTGCTGGACGGCTGCGGCTGCTGCCGCTCTGCGCCAAGCAGCTGGGAGAACTGTGCACGGAGCGTGAT
CCCTGCGACCCACACAAGGGTCTTCTGCGACTTCGGCTCCCCGCCAACCGAAGATTGGCGTGTGCA
CTGCCAAAGATGGTGCACCCTGTGTCTTCGGTGGGTCCGTGTACCGCAGCGGGCAGTCTTCCAAAGCAG
TTGCAATACCAAGTGCACCTGCCTGGATGGGGCCGTGGGCTGTGTGCCCTGTGCAGCATGGACGTGCGC
CTGCCCAGCCCTGACTGCCCTTCCCGAGAAGGGTCAAGCTGCCCGGAAATGCTGTGAGGAGTGGGTGT
GTGATGAGCCCAAGGACCGCACAGTGGTTGGCCCTGCCCTAGCTGCCTACCGACTGGAAGACACATTTGG
CCCTGACCCAATATGATGCGAGCCAATGCCTGGTCCAGACCACAGAGTGGAGCGCCTGTTCTAAGACC
TGTGGGATGGGCATCTCCACCCGGTTACCAATGACAATACCTTCTGCAGGCTGGAGAAGCAGAGTCGTC
TCTGCATGGTCAGGCCCTGTGAAGCTGACCTAGAGGAAAACATTAAGAAGGGCAAAAAGTGCATCCGGAC
GCCTAAAATTGCCAAGCCTGTCAAGTTTGAGCTTCTGGCTGCACCAAGTGTGAAGACCTACCGGGCTAAG
TTCTGTGGGTGTGCACGGACGGCCGCTGCTGCACCCGACAGAACACCACACTGCCGTTGGAGTTCA
AGTGCCCCGATGGCGAGATCATGAAAAAGAACATGATGTTTCATCAAGACCTGTGCCTGCCATTACAAGT
TCCCGGGGACAATGACATCTTTGAGTCTTGTACTACAGGAAGATGTATGGAGACATGGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR204597 representing NM_022266
Red=Cloning site Green=Tags(s)

MLASVAGPVSLALVLLLCTRPATGQDCSAQCQCAEAAPRCPAGVSLVLDGCGCCRVCAKQLGELCTERD
 PCDPHKGLFCDFGSPANRKIGVCTAKDGAPCVFGGSVYRSGESFQSSCKYQCTCLDGAVGCVPLCSMDVR
 LPSPDCPFPRRVKLPKCCCEEWVCDEPKDRTVVGPALAAYRLEDTFGPDPTMMRANCLVQTTEWSACSKT
 CGMGISTRVTNDNTFCRLEKQSRLCMVRPCEADLEENIKKGGKCI RTPKIAKPVKFELSGCTS VKTYRAK
 FCGVCTDGRCCTPHRTTLLPVEFKCPDGEIMKKNMMFIKTCACHYNC PGDNDIFESLYRKM YGDMA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_022266

ORF Size: 1041 bp

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_022266.2 , NP_071602.1
RefSeq Size:	2349 bp
RefSeq ORF:	1044 bp
Locus ID:	64032
UniProt ID:	Q9R1E9
Cytogenetics:	1p12
MW:	37.8 kDa
Gene Summary:	growth factor that may be a mediator of TGFbeta1-driven matrix production and normal skeletal remodeling [RGD, Feb 2006]