

Product datasheet for **RR204597L3V**

Ctgf (NM_022266) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ctgf (NM_022266) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ctgf
Synonyms:	CTGRP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_022266
ORF Size:	1041 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204597).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022266.2 , NP_071602.1
RefSeq Size:	2349 bp
RefSeq ORF:	1044 bp
Locus ID:	64032
UniProt ID:	Q9R1E9
Cytogenetics:	1p12
Gene Summary:	growth factor that may be a mediator of TGFbeta1-driven matrix production and normal skeletal remodeling [RGD, Feb 2006]



[View online »](#)