

Product datasheet for **RR204623L4V**

Sdc4 (NM_012649) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Sdc4 (NM_012649) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Sdc4
Synonyms:	RATRYUDOCA; RYUDOCA
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_012649
ORF Size:	606 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204623).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012649.2 , NP_036781.1
RefSeq Size:	2488 bp
RefSeq ORF:	609 bp
Locus ID:	24771
UniProt ID:	P34901
Cytogenetics:	3q42
Gene Summary:	regulates cell-matrix interactions; may bind phospholipids and play a role in phosphoinositide-mediated signaling [RGD, Feb 2006]



[View online »](#)