

Product datasheet for **RR204644L4V**

Adm2 (NM_201426) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Adm2 (NM_201426) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Adm2
Synonyms:	lmd
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_201426
ORF Size:	438 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204644).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_201426.1 , NP_958829.1
RefSeq Size:	1042 bp
RefSeq ORF:	441 bp
Locus ID:	399475
UniProt ID:	P61312
Cytogenetics:	7q34
Gene Summary:	acts via calcitonin receptor-like receptor (CRLR) and receptor activity-modifying protein (RAMP) protein complexes to induce cAMP production; decreases blood pressure; suppresses gastric emptying activity and food intake [RGD, Feb 2006]



[View online »](#)