

Product datasheet for **RR204667**

Glb1l (NM_001127529) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Glb1l (NM_001127529) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Glb1l
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>RR204667 representing NM_001127529
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCCACCCGACCTGCCTCCCTGCTCCTGCGGCTTGATTCTGCTGCTGCTGCCAGGCAGAAGCTC
GCTCCTTCGTAGTGGATCGGCAACATGACAGATTCCTCCTGGACGGGTCCCGTTCCGCTATATTTCTGG
CAGCCTGCACTACTTTCGCGTACCGCGGTTCTTTGGGCCGATCGACTTCTCAAGATGAAATTGAGCGGC
CTCAACGCTGTGCAATTTATGTGCCCTGGAACCTACCATGAGCCAGAACCTGGGATCTATAACTTTAATG
GCAGCCGGGACCTCACTGCCTTCTGAATGAGGCAGCTAAAGCAAACCTGTTGGTCATTCTGAGACCAGG
ACCTTACATCTGTGCTGAGTGGGAGATGGGAGGTCTCCCATCTTGGTTGCTTCGAAACCGTGATATTCAT
CTGAGAACCTCAGATCCTGCCTTCTTCATGCTGTGGACTCCTGGTTCAAGGTTTGTCTGCCAAAGATAT
ATCCATTTCTCTACCACAATGGGGCAACATCATTAGCATTGAGGTGGAGAATGAATATGGTAGCTACAA
GGTCTGTGATTTCAAGTATATGAGACACTTGGCTGGCCTCTCCGTGCATTACTAGGAGACGAGATCTTG
CTCTTCACCACAGATGGGCCTCAAGGACTCAGATGTGGTTCCTTCAGGGACTCTATGCCACTGTAGATT
TTGGCCCAGCTGACAACGTGACCAGAATCTTTCCCTGCTTCGGGAGTATGAACCTCATGGGCCATTGGT
GAACTCCGAGTATTATACAGGCTGGCTGGATTACTGGGGCAAAACCACTCCACACGGTCCACTCCAGCT
GTAGCTAGAGGACTAGAAAGAATACTCAAGTTGGGAGCCAGTGTGAACATGTATATGTTTCATGGAGGTA
CCAACCTTTGGATACTGGAATGGTCTGATGACAAGGGTACTTCCCTCCGATTACTACCAGCTACGACTA
TGATGCACCCATATCTGAAGCAGGAGACCCACACCTAAGCTTTTTCGCAATTCAAAATGCATCAGCAAG
TTCCAGGAAATTCCTTGGGACCTTACCTCCCCCAGCCCCAAGATGAAGCTTGGACCACTAACTATGA
GCCTGGATGGGAATTTGCTGCTTTCTAGACTTTCTCTGCCAGAAGGCCCATCCATTGATTTTGGC
ATTGACCTTTGAGGCTGTCAAACCTGGACCGTGGCTTTATGCTGTACCGGACCTATATGACTTACACTGTT
TCAGAACCAACACCATTGTGGTGCCAAATAATGGAATTCATGACCGTGCCTATGTGATGGTAGATGGGG
TATTGAAGGGCATTCTGGAACGGAGCTTGAAGAAGAAGTGTATTTGACTGGGAAAGTAGGGACCAGGCT
GGATATATTGCTAGAGAATATGGGGAGGCTCAGTTTTGGGTCTAACCACAGTGACTTCAAGGGCCTGTTA
GAACCACCACTTCTGGGGCAGACAGTCTCACTGGGTGGATGATGTTCCCTCTCAAGGTGGATAAGCTTG
TAAAATGGTGGATCCCCTTGCAGCTGATGAAAAGAGCCCAGCCTCAAGCTCCCTCTGCTCCCTGCCTTCTA
TTCTGCAAAGTTCCGGCTTAGGCTTACTTGGGGACACATTTCTCTATCTACCTGGATGGACCAAGGGC
CAAGTATGGATCAATGGGTTAACCTGGTCCGTACTGGACAATGCGGGTCCACAACAGACCCTATACG
TGCCAAGACCTCTGCTGTTTGGTAGCTCAATCAACAAAATCACATTGCTGGAACTAGAAGATGTGCCTCA
TGACCCCCAGATCCAATTTTGGATAAGCCTATCTCAATAGCACTTTGCATTGGGGGTATAACTATCTT
CTCTCAGAACCACAAGGTTCTATGAACCAATAGAGTTAAGTGGGCAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR204667 representing NM_001127529
 Red=Cloning site Green=Tags(s)


MPPDLPSLLLRLVILLLLQAARSFVVDQRHDFLLDGVPPFRYISGSLHYFRVPRVLWADRLKMKLSG
 LNAVQFYVPWNYHEPEPGIYNFNGSRDLTAFLNEAAKANLLVILRPGPYICAEWEMGLPSWLLRNRIH
 LRTSDPAFLHAVDSWFKVLLPKIYPFLYHNGGNIISIQVENEYGSYKVCDFKYMRLAGLFRALLGDEIL
 LFTTDGPQGLRCGSLQGLYATVDFGPADNVTRIFSLREYEPHGPLVNSEYTTGWLDYWGQNHSTRSTPA
 VARGLERILKLGASVNMVMFHGGTNFGYWNAGDDKGYFLPITTSYDYDAPISEAGDPTPKLFAIQNVISK
 FQEIPLGPLPPSPKMKLGPLTMSLDGNLLSFLDFLCPEGPIHSVLP LTFEAVKLDRGFMLYRTYMTYTV
 SEPTPLWVPNNGIHDRAYVMVDGVLKGI LERSLKEELYLTGKVGTRLDILLENMGRLSFGSNHSDFKGLL
 EPPLLGGTVLTGWMMFPLKVDKLVKWWIPLQLMKRAQPQAPSVP AFYSAKFPALGLLGDTFLYLPGWTKG
 QVWINGFNLGRYWTMRGPQQLYVPRPLLFSSINKITLLELEDVPHDPQIQFLDKPILNSTLHWGYNYL
 LSEPQGSYEPIELSGH


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:


ACCN: NM_001127529

ORF Size: 1938 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001127529.2](#), [NP_001121001.1](#)
RefSeq Size: 2334 bp
RefSeq ORF: 1941 bp
Locus ID: 301525
Cytogenetics: 9q33
MW: 73.3 kDa